

## RACK WASHER ACCESSORIES



Water Connection

Floor mounted quick-lock pipe





## Scientek General Purpose Wash Cart Model UC-6000

Scientek's general purpose wash cart is designed to hold all common sizes of rodent cages, debris pans and cage floors for thorough, efficient cleaning. This cart is built with heavy duty stainless steel components to withstand alkaline and acid detergents. Swivel-type retainers hold cages and pans firmly in place during the processing cycle.

## Scientek Bottle Washing Cart Model QLCT-6

Scientek's bottle washing cart connects easily to the optional floor mounted "quick-lock" fitting without needing to be tightened by hand. This makes loading our cart safe and simple. With capacity for up to six bottle baskets per load plus additional space for baskets of feeder bottle tops and spouts, this unit is a great help with productivity. A two position selector switch at the washer's control panel allows you to choose a special bottle washing spray pattern that ensures effective coverage from all angles.



## Scientek High Capacity Cage Washing System Model SMR-63



Scientek's equipment designers have created a high capacity double rack system for cleaning all common sizes of rodent cages. Our SMR-63 system holds up to 180 mouse cages or 60 rat cages in one load. This more than doubles the productivity of a standard 90" long rack washer. Along with the two cage holding racks we furnish a central spray header to distribute the wash and rinse solutions to all parts of the load. The central spray header connects easily to the optional floor mounted "quick-lock" fitting without needing to be tightened by hand. All components are stainless steel and acid compatible.

## **Scientek Spindle Type Bottle Washing Cart Model QLSP-6**



Scientek's spindle-type bottle washing cart connects easily to the optional floor mounted "quick-lock" fitting without needing to be tightened by hand. This makes loading our carts safe and simple. This speciality cart cleans up to 108 feeder bottles, tops and spouts in one load. A two position selector switch at the washers control panel allows you to choose a special bottle washing spray pattern to ensure effective coverage from all angles, including the interior of each bottle. All components are stainless steel and acid compatible.

**Scientek Cage & Rack Washer Accessory Dimensions & Weights**

<b>Model #</b>	<b>SAE Dimensions</b>	<b>Metric Dimensions</b>	<b>Weight</b>
SMR-63	43.5" W x 86" H x 84" L	1105 mm x 2184 mm x 2134 mm	625 Lbs - 1375 Kg
UC-6000	32" W x 67" H x 77.5" L	813 mm x 1702 mm x 1969 mm	452 Lbs - 995 Kg
QLCT-6	39" W x 42" H x 64" L	991 mm x 1067 mm x 1626 mm	120 Lbs - 264 Kg
QLSP-6	39" W x 42" H x 64" L	991 mm x 1067 mm x 1626 mm	200 Lbs - 440 Kg

**Note:**

Machines using the accessories SMR-63, QLCT-6 or QLSP-6 must have the Scientek quick-lock connection in the washer floor and may require an additional stainless steel guide track in the cabinet ceiling.

Machines using the SMR-63 must have a minimum cabinet interior size of 48" wide x 88" high x 90" long.